


Complex PCR set-up simplified.

APPLICATION PACKAGE - FREEDOM EVO 75

Effortless automated PCR set-up on the Freedom EVO® 75 workstation, with a simple touchscreen-based operator interface.

Fully automated reaction set-up for complex PCR plate layouts, with up to 96 master mixes and multiple PCR plates per run, plus individual controls and dilution curves. A simple software wizard guides you through instrument set-up, including defining and reviewing your plate layout and starting the run.

Get the full application package, including the liquid handling platform, software wizard and training.



What can I use it for?


The Freedom EVO 75 PCR workstation and TouchTools™ PCR Wizard simplifies reaction set-up for a wide range of applications – from end-point, real-time and multiplex PCR protocols, to genotyping and gene expression methods, as well as pathogen and mutagenesis detection.

Flexible and user-friendly five-step PCR set-up:

1. Select samples
2. Verify PCR plate layout
3. Load worktable
4. Watch progress information
5. Export data to qPCR instrument

Intuitive TouchTools interface:

- Visual prompts for loading reagents, samples and labware, including required quantities/volumes
- Helps to reduce training time and costs




Step-by-step
guided
instrument
set-up

Seamless
data transfer
to real-time
cyclers

The Freedom EVO PCR workstation includes:

- 1 Freedom EVO 75
- 2 Liquid handling arm
- 3 Storage for disposable tips
- 4 Waste chute for disposable tips
- 5 Tube racks
- 6 Plate carriers
- 7 Trough racks
- 8 TouchTools™ PCR Wizard software
 - Computer with touchscreen monitor
 - Freedom EVOware® Standard software



www.tecan.com

 **TECAN.**

Australia +61 3 9647 4100 **Austria** +43 62 46 89 33 **Belgium** +32 15 42 13 19 **China** +86 21 220 63 206 **Denmark** +45 70 23 44 50 **France** +33 4 72 76 04 80
Germany +49 79 51 94 170 **Italy** +39 02 92 44 790 **Japan** +81 44 556 73 11 **Netherlands** +31 18 34 48 17 4 **Singapore** +65 644 41 886 **Spain** +34 93 490 01 74
Sweden +46 8 750 39 40 **Switzerland** +41 44 922 89 22 **UK** +44 118 9300 300 **USA** +1 919 361 5200 **Other countries** +43 62 46 89 33

© 2016 Tecan trading AG, Switzerland, all rights reserved. For disclaimer and trademarks please visit www.tecan.com

